
Product Name: ZHX2 Rabbit Monoclonal Antibody**Catalog #: AMRe85168**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200
Molecular Weight	Calculated MW: 92 kDa; Observed MW: 110 kDa

Antigen Information


Gene Name	ZHX2
Alternative Names	AFP regulator 1; Zinc finger and homeodomain protein 2
Gene ID	22882.0
SwissProt ID	Q9Y6X8
Immunogen	Recombinant protein of human ZHX2

Background


Acts as a transcriptional repressor.

Research Area

Image Data



Western blot analysis of ZHX2 in rat Brain, C6, HeLa lysates using ZHX2 antibody.



Immunocytochemistry analysis of ZHX2 (green) in MCF-7 using ZHX2 antibody, and DAPI (blue)