

Product Name: StAR Rabbit Monoclonal Antibody**Catalog #: AMRe85190**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 32 kDa

Antigen Information


Gene Name	StAR
Alternative Names	STARD1
Gene ID	6770.0
SwissProt ID	P49675
Immunogen	A synthetic peptide of human StAR

Background

Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer membrane to the inner mitochondrial membrane.

Research Area

Image Data



Western blot analysis of StAR in K562, 3T3, HeLa lysates using StAR antibody.