

---

**Product Name: R Cadherin Rabbit Monoclonal Antibody****Catalog #: AMRe85212**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IP 1:10-1:20
<b>Molecular Weight</b>	Calculated MW: 100 kDa; Observed MW: 130 kDa

**Antigen Information**


<b>Gene Name</b>	R Cadherin
<b>Alternative Names</b>	cadherin 4; type 1; R-cadherin (retinal); CAD4; RCAD; R-CAD
<b>Gene ID</b>	1002.0
<b>SwissProt ID</b>	P55283
<b>Immunogen</b>	A synthetic peptide of human R Cadherin

**Background**


Cadherins are calcium-dependent cell adhesion proteins.

**Research Area**

## Image Data



Western blot analysis of RCadherin/Cadherin 4 in mouse brain, U87-MG lysates using R Cadherin antibody.



Western blot analysis of RCadherin/Cadherin 4 in rat Brain, C6 lysates using RCadherin/Cadherin 4 antibody.