

Product Name: ACADS Rabbit Monoclonal Antibody**Catalog #: AMRe85247**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Antigen Information


Gene Name	ACADS
Alternative Names	SCAD; ACAD3
Gene ID	35.0
SwissProt ID	P16219
Immunogen	A synthetic peptide of human ACADS

Background


Introduces a double bond at position 2 in saturated acyl-CoA's of short chain length, i.e. less than 6 carbon atoms.

Research Area

Image Data



Western blot analysis of ACADS in mouse liver, 3T3 lysates using ACADS antibody.



Western blot analysis of ACADS in Jurkat, C6, HeLa lysates using ACADS antibody.