
Product Name: AER61 Rabbit Monoclonal Antibody**Catalog #: AMRe85263**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Antigen Information


Gene Name	AER61
Alternative Names	Extracellular O-linked N-acetylglucosamine transferase
Gene ID	285203.0
SwissProt ID	Q5NDL2
Immunogen	Recombinant protein of human AER61

Background

EGF domain-specific O-linked N-acetylglucosamine transferase

Research Area

Image Data



Western blot analysis of AER61 in C6, HeLa lysates using AER61 antibody.