
Product Name: ALDH1L1 Rabbit Monoclonal Antibody**Catalog #: AMRe85269**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 99 kDa; Observed MW: 99 kDa

Antigen Information


Gene Name	ALDH1L1
Alternative Names	FDH; FTHFD; 10-fTHF; 10-FTHFDH
Gene ID	10840.0
SwissProt ID	O75891
Immunogen	A synthetic peptide of human ALDH1L1

Background

Highly expressed in liver, pancreas and kidney.

Research Area

Image Data



Western blot analysis of ALDH1L1 in rat Brain, 3T3 lysates using ALDH1L1 antibody.