
Product Name: Angiotensinogen Rabbit Monoclonal Antibody**Catalog #: AMRe85286**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa

Antigen Information


Gene Name	Angiotensinogen
Alternative Names	Angiotensinogen; Serpin A8; AGT; SERPINA8; Des-Asp[1]-angiotensin II; Serine (or cysteine) proteinase inhibitor
Gene ID	11606.0
SwissProt ID	P11859
Immunogen	Recombinant protein of mouse Angiotensinogen

Background

Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis.

Research Area

Image Data



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.