

Product Name: Annexin VII Rabbit Monoclonal Antibody
Catalog #: AMRe85289

Summary

Production Name	Annexin VII Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,IHC
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Annexin VII
Alternative Names	SNX; ANX7; SYNEXIN
Gene ID	310.0
SwissProt ID	P20073.Recombinant protein of human Annexin VII

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:10-1:20,IHC:1:50-1:100
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 51,47 kDa

Background


Product Name: Annexin VII Rabbit Monoclonal Antibody
Catalog #: AMRe85289




Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Research Area

Image Data



Western blot analysis of Annexin VII in K562, rat Brain, HeLa lysates using Annexin VII antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Annexin VII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.