
Product Name: Argonaute 2 Rabbit Monoclonal Antibody**Catalog #: AMRe85304**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200,IP 1:10-1:20
Molecular Weight	Calculated MW: 97 kDa; Observed MW: 97 kDa

Antigen Information

Gene Name	Argonaute 2
Alternative Names	Ago 2; Argonaute 2; dAgo2; eIF2C 2; hAgo2; PPD
Gene ID	27161.0
SwissProt ID	Q9UKV8
Immunogen	A synthetic peptide of human Argonaute-2


Background

Required for RNA-mediated gene silencing (RNAi) by the RNA-induced silencing complex (RISC). Can also upregulate the translation of specific mRNAs under certain growth conditions. Binds to the AU element of the 3'-UTR of the TNF (TNF-alpha) mRNA and upregulates translation under conditions of serum starvation. Also required for transcriptional gene silencing (TGS),


in which short RNAs known as antigen RNAs or agRNAs direct the transcriptional repression of complementary promoter regions.

Research Area

Image Data



Western blot analysis of Argonaute2 in K562, MCF-7, HeLa lysates using Argonaute2 antibody.



Immunocytochemistry analysis of Argonaute 2 (green) in HCT116 using Argonaute 2 antibody, and DAPI (blue).