

Product Name: ARHGAP25 Rabbit Monoclonal Antibody**Catalog #: AMRe85305**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Antigen Information


Gene Name	ARHGAP25
Alternative Names	Rho GTPase activating protein 25; KAIA0053
Gene ID	9938.0
SwissProt ID	P42331
Immunogen	Recombinant protein of human ARHGAP25

Background

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Research Area

Image Data



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.