
Product Name: ATP1B3 Rabbit Monoclonal Antibody**Catalog #: AMRe85325**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 32 kDa; Observed MW: 40 kDa

Antigen Information


Gene Name	ATP1B3
Alternative Names	atp1b3; ATPB3; CD298; FLJ29027; NKAB3S; Sodium pump subunit beta 3; Sodium/potassium dependent ATPase subunit beta 3
Gene ID	483.0
SwissProt ID	P54709
Immunogen	A synthetic peptide of human ATP1B3

Background

This is the non-catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of Na(+) and K(+) ions across the plasma membrane. The exact function of the beta-3 subunit is not known.

Research Area

Image Data



Western blot analysis of ATP1B3 in C6, 3T3 lysates using ATP1B3 antibody.