
Product Name: AVEN Rabbit Monoclonal Antibody**Catalog #: AMRe85330**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 50 kDa

Antigen Information


Gene Name	AVEN
Alternative Names	PDCD12
Gene ID	57099.0
SwissProt ID	Q9NQS1
Immunogen	A synthetic peptide of human AVEN

Background


Protects against apoptosis mediated by Apaf-1.

Research Area

Image Data



Western blot analysis of AVEN in rat Brain, MCF-7, 3T3, HeLa lysates using AVEN antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using AVEN antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.