

Product Name: BAP31 Rabbit Monoclonal Antibody
Catalog #: AMRe85341



Summary

Production Name	BAP31 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	BAP31
Alternative Names	BCAP31; BAP31; DXS1357E; B-cell receptor-associated protein 31; BCR-associated protein 31; Bap31; 6C6-AG tumor-associated antigen; Protein CDM; p28
Gene ID	10134.0
SwissProt ID	P51572.A synthetic peptide of human BAP31

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:100,ICC:1:50-1:200
Molecular Weight	Calculated MW: 28 kDa; Observed MW: 28 kDa

Product Name: BAP31 Rabbit Monoclonal Antibody
Catalog #: AMRe85341




Background

May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi. May be involved in CASP8-mediated apoptosis.


Research Area

Apoptosis


Image Data



Western blot analysis of BAP31 in K562, HeLa lysates using BAP31 antibody.



Immunocytochemistry analysis of BAP31(green) in HeLa using BAP31 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using BAP31 antibody. High-pressure and

Product Name: BAP31 Rabbit Monoclonal Antibody
Catalog #: AMRe85341



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.