

Product Name: Bcl10 Rabbit Monoclonal Antibody
Catalog #: AMRe85344



Summary

Production Name	Bcl10 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Bcl10 BCL10; CIPER; CLAP; B-cell lymphoma/leukemia 10; B-cell CLL/lymphoma 10; Bcl-10;
Alternative Names	CARD-containing molecule enhancing NF-kappa-B; CARD-like apoptotic protein; hCLAP; CED-3/ICH-1 prodomain homologous E10-like regulator; CIPER; Cellular homolog
Gene ID	8915.0
SwissProt ID	O95999. A synthetic peptide of human Bcl10

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 26 kDa; Observed MW: 26 kDa

Product Name: Bcl10 Rabbit Monoclonal Antibody
Catalog #: AMRe85344




Background

Promotes apoptosis, pro-caspase-9 maturation and activation of NF-kappa-B via NIK and IKK. May be an adapter protein between upstream TNFR1-TRADD-RIP complex and the downstream NIK-IKK-IKAP complex. Is a substrate for MALT1.

Research Area

Apoptosis

Image Data



Western blot analysis of Bcl10 in K562, C6, 3T3, HeLa lysates using Bcl10 antibody.

Note

For research use only.