

Product Name: BRD2 Rabbit Monoclonal Antibody
Catalog #: AMRe85360



Summary

Production Name	BRD2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	BRD2
Alternative Names	Brd2; FSH; FSRG1; NAT; O27.1.; RING3; RNF3
Gene ID	6046.0
SwissProt ID	P25440. A synthetic peptide of human BRD2

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 110 kDa

Background


May play a role in spermatogenesis or folliculogenesis.

Product Name: BRD2 Rabbit Monoclonal Antibody
Catalog #: AMRe85360



Research Area

Image Data



Western blot analysis of BRD2 in rat Brain, Jurkat lysates using BRD2 antibody.

Note

For research use only.