

Product Name: GC1q R Rabbit Monoclonal Antibody**Catalog #: AMRe85364**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 31 kDa; Observed MW: 31 kDa

Antigen Information

Gene Name	GC1q R C1QBP; GC1QBP; HABP1; SF2P32; Complement component 1 Q subcomponent-binding
Alternative Names	protein; mitochondrial; GC1q-R protein; Glycoprotein gC1qBP; C1qBP; Hyaluronan-binding protein 1; Mitochondrial matrix protein p32; p33
Gene ID	708.0
SwissProt ID	Q07021
Immunogen	A synthetic peptide of human C1QBP


Background

Binds to the globular "heads" of C1Q thus inhibiting C1 activation.

Research Area

Apoptosis

Image Data



Western blot analysis of GC1q R in K562, rat Brain, C6, 3T3, HeLa lysates using GC1q R antibody.