
Product Name: C3orf38 Rabbit Monoclonal Antibody**Catalog #: AMRe85366**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Antigen Information


Gene Name	C3orf38
Alternative Names	Uncharacterized protein C3orf38
Gene ID	285237.0
SwissProt ID	Q5JPI3
Immunogen	Recombinant protein of human C3orf38

Background

May be involved in apoptosis regulation.

Research Area

Image Data



Western blot analysis of C3orf38 in K562, HeLa lysates using C3orf38 antibody.