
Product Name: CARD9 Rabbit Monoclonal Antibody**Catalog #: AMRe85377**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Antigen Information


Gene Name	CARD9
Alternative Names	Caspase recruitment domain-containing protein 9; hCARD9
Gene ID	64170.0
SwissProt ID	Q9H257
Immunogen	A synthetic peptide of human CARD9

Background

Caspase recruitment domain-containing protein 9

Research Area

Image Data



Western blot analysis of CARD9 in K562, THP-1 lysates using CARD9 antibody.