

Product Name: CD134/OX40L Receptor Rabbit Monoclonal Antibody**Catalog #: AMRe85394**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 45 kDa

Antigen Information


Gene Name	CD134/OX40L Receptor
Alternative Names	OX40 antigen; OX40L receptor
Gene ID	22163.0
SwissProt ID	P47741
Immunogen	Recombinant protein of mouse CD134

Background

Receptor for TNFSF4/OX40L/GP34. Is a costimulatory molecule implicated in long-term T-cell immunity (By similarity).

Research Area

Image Data



Western blot analysis of CD134 in mouse spleen lysates using CD134/OX40L Receptor antibody.