

Product Name: CD35 Rabbit Monoclonal Antibody
Catalog #: AMRe85406



Summary

Production Name	CD35 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	CD35
Alternative Names	C3 binding protein; C3b/C4b receptor; C3BR; C4BR; CD35; CR1; KN
Gene ID	1378.0
SwissProt ID	P17927. A synthetic peptide of human CD35

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:100
Molecular Weight	Calculated MW: 224 kDa; Observed MW: 280 kDa

Background


Mediates cellular binding of particles and immune complexes that have activated complement.

Product Name: CD35 Rabbit Monoclonal Antibody
Catalog #: AMRe85406



Research Area

Image Data



Western blot analysis of CD35 in Human tonsil lysates using CD35 antibody.

Note

For research use only.