
Product Name: CD3G Rabbit Monoclonal Antibody**Catalog #: AMRe85409**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Antigen Information


Gene Name	CD3G
Alternative Names	CD3 gamma; CD3g; T3G
Gene ID	917.0
SwissProt ID	P09693
Immunogen	A synthetic peptide of human CD3G

Background


The CD3 complex mediates signal transduction.

Research Area

Image Data



Western blot analysis of CD3G in Jurkat lysates using CD3G antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD3G antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.