

**Product Name: Cellular Apoptosis Susceptibility Rabbit
Monoclonal Antibody
Catalog #: AMRe85427**



Summary

Production Name	Cellular Apoptosis Susceptibility Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Cellular Apoptosis Susceptibility
Alternative Names	CSE1L; CAS; CSE1; Exp2; Importin-alpha re-exporter; Exportin-2; XPO2
Gene ID	1434.0
SwissProt ID	P55060. A synthetic peptide of human Cellular Apoptosis Susceptibility

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:100,ICC:1:50-1:200
Molecular Weight	Calculated MW: 110 kDa; Observed MW: 110 kDa

Background


**Product Name: Cellular Apoptosis Susceptibility Rabbit
Monoclonal Antibody
Catalog #: AMRe85427**



Export receptor for importin-alpha. Mediates importin-alpha re-export from the nucleus to the cytoplasm after import substrates (cargos) have been released into the nucleoplasm. In the nucleus binds cooperatively to importin-alpha and to the GTPase Ran in its active GTP-bound form.

Research Area

Image Data



Western blot analysis of ular Apoptosis Susceptibility in 3T3, HeLa lysates using ular Apoptosis Susceptibility antibody.

Note

For research use only.