

Summary

Production Name	ACAP2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	ACAP2
Alternative Names	CENTB2; CNT-B2
Gene ID	23527.0
SwissProt ID	Q15057. A synthetic peptide of human Centaurin beta 2

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 88 kDa

Background


GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).

Product Name: ACAP2 Rabbit Monoclonal Antibody
Catalog #: AMRe85429



Research Area

Image Data



Western blot analysis of ACAP2 in HeLa, CHO-K1, C6, Jurkat lysates using ACAP2 antibody.

Note

For research use only.