

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe85434



Summary

Production Name	Chromogranin A Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Chromogranin A
Alternative Names	beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I; Vasostatin-2; Vasostatin II; Pancreastatin; CHGA
Gene ID	1113.0
SwissProt ID	P10645. A synthetic peptide of human Chromogranin A

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:100,ICC:1:50-1:200
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 92 kDa

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe85434




Background

Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

Research Area

Image Data



Western blot analysis of Chromogranin A in PC-12 lysates using Chromogranin A antibody.

Note

For research use only.