

Product Name: CLPTM1 Rabbit Monoclonal Antibody
Catalog #: AMRe85441



Summary

Production Name	CLPTM1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	CLPTM1
Alternative Names	Cleft lip and palate transmembrane protein 1
Gene ID	1209.0
SwissProt ID	O96005. A synthetic peptide of human CLPTM1

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 76 kDa; Observed MW: 90 kDa

Background


May play a role in T-cell development.

Product Name: CLPTM1 Rabbit Monoclonal Antibody
Catalog #: AMRe85441



Research Area

Image Data



Western blot analysis of CLPTM1 in 293T, U87-MG, Jurkat lysates using CLPTM1 antibody.

Note

For research use only.