
Product Name: COG6 Rabbit Monoclonal Antibody**Catalog #: AMRe85443**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Antigen Information


Gene Name	COG6
Alternative Names	COD2; SHNS; CDG2L
Gene ID	57511.0
SwissProt ID	Q9Y2V7
Immunogen	A synthetic peptide of human COG6

Background

Required for normal Golgi function.

Research Area

Image Data



Western blot analysis of COG6 in K562, 3T3 lysates using COG6 antibody.