
Product Name: CSHL1 Rabbit Monoclonal Antibody**Catalog #: AMRe85457**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Antigen Information


Gene Name	CSHL1
Alternative Names	Lactogen-like; Chorionic somatomammotropin-like
Gene ID	1444.0
SwissProt ID	Q14406
Immunogen	Recombinant protein of human CSHL1

Background

May be a novel gestational hormone required to compensate for absence of other members of the GH/CS cluster during gestation.

Research Area

Image Data



Western blot analysis of CSHL1 in Human placenta lysates using CSHL1 antibody.