
Product Name: DAP1 Rabbit Monoclonal Antibody**Catalog #: AMRe85487**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

Antigen Information

Gene Name	DAP1
Alternative Names	Death-associated protein 1; DAP-1
Gene ID	1611.0
SwissProt ID	P51397
Immunogen	A synthetic peptide of human DAP1


Background

Negative regulator of autophagy. Involved in mediating interferon-gamma-induced cell death.

Research Area

Autophagy, Apoptosis

Image Data



Western blot analysis of DAP1 in MCF-7, C6, 3T3 lysates using DAP1 antibody.