
Product Name: DEK Rabbit Monoclonal Antibody**Catalog #: AMRe85500**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 50 kDa

Antigen Information


Gene Name	DEK
Alternative Names	Dek; DEK oncogene; Protein DEK
Gene ID	7913.0
SwissProt ID	P35659
Immunogen	Recombinant protein of human DEK

Background


Involved in chromatin organization.

Research Area

Image Data



Western blot analysis of DEK in K562, MCF-7, HeLa lysates using DEK antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DEK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.