

Product Name: GSDME Rabbit Monoclonal Antibody
Catalog #: AMRe85502



Summary

Production Name	GSDME Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	GSDME
Alternative Names	DFNA5; ICERE-1; GSDME
Gene ID	1687.0
SwissProt ID	O60443. Recombinant protein of human DFNA5/GSDME

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa

Background

Plays a role in the TP53-regulated cellular response to DNA damage probably by cooperating with TP53


Product Name: GSDME Rabbit Monoclonal Antibody
Catalog #: AMRe85502



(PubMed:16897187, PubMed:18223688).

Research Area

Image Data



Western blot analysis of DFNA5/GSDME in A549 lysates using GSDME antibody.

Note

For research use only.