

Product Name: DHX8 Rabbit Monoclonal Antibody**Catalog #: AMRe85506**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC 1:50-1:200
Molecular Weight	Calculated MW: 139 kDa; Observed MW: 139 kDa

Antigen Information


Gene Name	DHX8
Alternative Names	DDX8; Dhr2; HRH1; PRP22; PRPF22
Gene ID	1659.0
SwissProt ID	Q14562
Immunogen	A synthetic peptide of human DHX8

Background

Facilitates nuclear export of spliced mRNA by releasing the RNA from the spliceosome.

Research Area

Image Data



Western blot analysis of DHX8 in HeLa, 293T lysates using DHX8 antibody.