
Product Name: DLX4 Rabbit Monoclonal Antibody**Catalog #: AMRe85510**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200,IP 1:10-1:20
Molecular Weight	Calculated MW: 26 kDa; Observed MW: 35 kDa

Antigen Information


Gene Name	DLX4
Alternative Names	DLX4; BP1; DLX7; DLX8; DLX9; Homeobox protein DLX-4; Beta protein 1; Homeobox protein DLX-7; Homeobox protein DLX-8
Gene ID	1748.0
SwissProt ID	Q92988
Immunogen	A synthetic peptide of human Dlx4

Background


May play a role in determining the production of hemoglobin S. May act as a repressor. During embryonic development, plays a role in palatogenesis.

Research Area

Image Data



Western blot analysis of Dlx4 in K562, rat Brain, HeLa lysates using Dlx4 antibody.



Immunocytochemistry analysis of DLX4 (green) in K562 using DLX4 antibody, and DAPI (blue).

