
Product Name: DNAJA4 Rabbit Monoclonal Antibody**Catalog #: AMRe85512**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Antigen Information


Gene Name	DNAJA4
Alternative Names	MST104; MSTP104; PRO1472
Gene ID	55466.0
SwissProt ID	Q8WW22
Immunogen	A synthetic peptide of human DNAJA4

Background


Enables chaperone binding activity and unfolded protein binding activity. Involved in several processes, including negative regulation of endothelial cell migration; negative regulation of inclusion body assembly; and protein refolding. Located in cytosol and membrane.

Research Area

Image Data



Western blot analysis of DNAJA4 in HeLa lysates using DNAJA4 antibody.



Western blot analysis of DNAJA4 in C6, rat Brain lysates using DNAJA4 antibody.