

Product Name: DPS1 Rabbit Monoclonal Antibody**Catalog #: AMRe85519**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa

Antigen Information


Gene Name	DPS1
Alternative Names	Trans-prenyltransferase 1; TPT 1; Solanesyl-diphosphate synthase subunit 1; Decaprenyl-diphosphate synthase subunit 1
Gene ID	23590.0
SwissProt ID	Q5T2R2
Immunogen	A synthetic peptide of human DPS1

Background

Decaprenyl-diphosphate synthase subunit 1

Research Area

Image Data



Western blot analysis of DPS1 in K562 lysates using DPS1 antibody.