

Product Name: EN1 Rabbit Monoclonal Antibody
Catalog #: AMRe85537



Summary

Production Name	EN1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	EN1
Alternative Names	EN1; Engrailed 1; Engrailed; Hu-En-1; Mo-en.1
Gene ID	2019.0
SwissProt ID	Q05925. A synthetic peptide of human EN1/Engrailed 1

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 50 kDa

Background


Belongs to the engrailed homeobox family.

Product Name: EN1 Rabbit Monoclonal Antibody
Catalog #: AMRe85537



Research Area

Image Data



Western blot analysis of EN1/Engrailed 1 in rat Brain lysates using EN1 antibody.

Note

For research use only.