

Product Name: ENO3 Rabbit Monoclonal Antibody**Catalog #: AMRe85538**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Antigen Information


Gene Name	ENO3
Alternative Names	enolase 3; MSE; GSD13
Gene ID	2027.0
SwissProt ID	P13929
Immunogen	A synthetic peptide of human ENO3

Background


Appears to have a function in striated muscle development and regeneration.

Research Area

Image Data



Western blot analysis of ENO3 in rat Brain, C6, 3T3, HeLa lysates using ENO3 antibody.



Immunocytochemistry analysis of ENO1/ENO2/ENO3(green) in HeLa using ENO1/ENO2/ENO3 antibody, and DAPI(blue)