

Product Name: Fetuin A Rabbit Monoclonal Antibody**Catalog #: AMRe85564**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 45-50 kDa

Antigen Information


Gene Name	Fetuin A
Alternative Names	16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743
Gene ID	197.0
SwissProt ID	P02765
Immunogen	A synthetic peptide of human Fetuin A

Background

Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Research Area

Image Data



Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.