

Product Name: FKBP51 Rabbit Monoclonal Antibody
Catalog #: AMRe85569



Summary

Production Name	FKBP51 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	FKBP51
Alternative Names	P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10
Gene ID	2289.0
SwissProt ID	Q13451. A synthetic peptide of human FKBP51

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 51 kDa

Background


Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Product Name: FKBP51 Rabbit Monoclonal Antibody
Catalog #: AMRe85569



Research Area

Image Data



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.

Note

For research use only.