

Product Name: Flotillin 1 Rabbit Monoclonal Antibody**Catalog #: AMRe85571**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Antigen Information


Gene Name	Flotillin 1
Alternative Names	Flotillin-1
Gene ID	10211.0
SwissProt ID	O75955
Immunogen	A synthetic peptide of human Flotillin 1

Background


May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Research Area

Image Data



Western blot analysis of Flotillin 1 in K562, HeLa lysates using Flotillin 1 antibody.



Immunocytochemistry analysis of Flotillin 1(green) in HeLa using Flotillin 1 antibody, and DAPI(blue)