

Product Name: FNDC5 Rabbit Monoclonal Antibody**Catalog #: AMRe85572**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 24 kDa; Observed MW: 21 kDa

Antigen Information


Gene Name	FNDC5
Alternative Names	FNDC 5; FRCP2; irisin
Gene ID	252995.0
SwissProt ID	Q8NAU1
Immunogen	A synthetic peptide of human FNDC5

Background

Irisin: Contrary to mouse, may not be involved in the beneficial effects of muscular exercise, nor in the induction of browning of human white adipose tissue.

Research Area

Image Data



Western blot analysis of FNDC5 in K562, C6, 3T3, HeLa lysates using FNDC5 antibody.