

Product Name: Follistatin Rabbit Monoclonal Antibody
Catalog #: AMRe85576



Summary

Production Name	Follistatin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Follistatin
Alternative Names	Activin binding protein; Follistatin; FS; fst; Human follistatin gene
Gene ID	10468.0
SwissProt ID	P19883. Recombinant protein of human Follistatin

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Background

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of


Product Name: Follistatin Rabbit Monoclonal Antibody
Catalog #: AMRe85576



pituitary follicle stimulating hormone (FSH).

Research Area

Image Data



Western blot analysis of Follistatin in rat Brain lysates using Follistatin antibody.

Note

For research use only.