

Product Name: FYB Rabbit Monoclonal Antibody**Catalog #: AMRe85588**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 85 kDa; Observed MW: 130,120 kDa

Antigen Information


Gene Name	FYB
Alternative Names	ADAP; Fyb; FYN binding protein; p120/p130; SLAP130
Gene ID	2533.0
SwissProt ID	O15117
Immunogen	A synthetic peptide of human FYB

Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

Research Area

Image Data



Western blot analysis of FYB in Jurkat lysates using FYB antibody.