

Product Name: GLI3 Rabbit Monoclonal Antibody**Catalog #: AMRe85613**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC 1:50-1:200
Molecular Weight	Calculated MW: 170 kDa; Observed MW: 200 kDa

Antigen Information


Gene Name	GLI3
Alternative Names	GLI3-190; GLI3FL; GLI3-83; GLI3 C-terminally truncated form
Gene ID	463369.0
SwissProt ID	Q5IS56
Immunogen	Recombinant protein of Gli3

Background

Has a dual function as a transcriptional activator and a repressor of the sonic hedgehog (Shh) pathway, and plays a role in limb development.

Research Area

Image Data



Western blot analysis of Gli3 in 293 lysates using GLI3 antibody.