

Product Name: Glutathione Reductase Rabbit Monoclonal Antibody**Catalog #: AMRe85618**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 51 kDa

Antigen Information


Gene Name	Glutathione Reductase
Alternative Names	GR; HEL-75; HEL-S-122m
Gene ID	2936.0
SwissProt ID	P00390
Immunogen	A synthetic peptide of human Glutathione Reductase

Background

Maintains high levels of reduced glutathione in the cytosol.

Research Area

Image Data



Western blot analysis of Glutathione Reductase in HeLa, Jurkat lysates using Glutathione Reductase antibody.