
Product Name: GSTK1 Rabbit Monoclonal Antibody**Catalog #: AMRe85619**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,ICC 1:50-1:200,IP 1:10-1:20
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Antigen Information


Gene Name	GSTK1
Alternative Names	GST; GST class kappa; GST13; Gstk1; hGSTK1
Gene ID	373156.0
SwissProt ID	Q9Y2Q3
Immunogen	A synthetic peptide of human Glutathione S Transferase kappa 1

Background


Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene (CDNB).

Research Area

Image Data



Western blot analysis of GSTK1 in HeLa lysates using GSTK1 antibody.



Immunocytochemistry analysis of Glutathione S Transferase kappa 1 (green) in HeLa using Glutathione S Transferase kappa 1 antibody, and DAPI (blue)