
Product Name: GPRC5B Rabbit Monoclonal Antibody**Catalog #: AMRe85632**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:10-1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Antigen Information


Gene Name	GPRC5B
Alternative Names	RAIG2; RAIG-2
Gene ID	51704.0
SwissProt ID	Q9NZH0
Immunogen	A synthetic peptide of human GPRC5B

Background


G-protein coupled receptor family C group 5 member B

Research Area

Image Data



Western blot analysis of GPRC5B in A549 lysates using GPRC5B antibody.



Western blot analysis of GPRC5B in rat Brain, 3T3 lysates using GPRC5B antibody.