

Product Name: HARS Rabbit Monoclonal Antibody**Catalog #: AMRe85637**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, IP 1:10-1:20
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 57 kDa

Antigen Information


Gene Name	HARS
Alternative Names	HRS; CMT2W; USH3B
Gene ID	3035.0
SwissProt ID	P12081
Immunogen	A synthetic peptide of human HARS

Background

Cytoplasmic histidine--tRNA ligase (Probable). Plays a role in axon guidance.

Research Area

Image Data



Western blot analysis of HARS in K562, HeLa lysates using HARS antibody.