
Product Name: HNRNPA1L2 Rabbit Monoclonal Antibody**Catalog #: AMRe85645**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Antigen Information


Gene Name	HNRNPA1L2
Alternative Names	hnRNP core protein A1-like 2
Gene ID	144983.0
SwissProt ID	Q32P51
Immunogen	A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)

Background

Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.

Research Area

Image Data



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.