

Product Name: HLA C Rabbit Monoclonal Antibody
Catalog #: AMRe85658



Summary

Production Name	HLA C Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide, 0.05% BSA and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	HLA C
Alternative Names	MHC; HLAC; HLC-C; D6S204; PSORS1; HLA-JY3
Gene ID	3107.0
SwissProt ID	P10321. A synthetic peptide of human HLA-C

Application

Dilution Ratio	WB:1:500-1:1000, IP:1:10-1:20
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Background


Antigen-presenting major histocompatibility complex class I (MHCI) molecule with an important role in reproduction and antiviral immunity.

Product Name: HLA C Rabbit Monoclonal Antibody
Catalog #: AMRe85658



Research Area

Image Data



Western blot analysis of HLAC in Hela, 293T, Jurkat, 3T3, HepG2 lysates using HLA C antibody.

Note

For research use only.