

Product Name: HNF 4 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe85663

Summary

Production Name	HNF 4 alpha Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide,0.05%BSA and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	HNF 4 alpha
Alternative Names	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha
Gene ID	15378.0
SwissProt ID	P49698. Recombinant protein of mouse HNF-4-alpha

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:100,ICC:1:50-1:200
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Background


Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin,

Product Name: HNF 4 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe85663

apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and inTine.MiscellaneousBinds fatty acids.

Research Area

Image Data



Western blot analysis of HNF4alpha in Raw264.7 lysates using HNF 4 alpha antibody.

Note

For research use only.